

STAT1 (6C2) Mouse mAb


ABZ3971

Product Overview


Name	STAT1 (6C2) Mouse mAb
Accession(Primary)	P42224
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 91 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P42224)

Gene Name	STAT1
About Gene Symbol	STAT1 (Signal Transducer And Activator Of Transcription 1) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cell



Western blot analysis of STAT1 in HeLa, MCF-7, COS7, SW480 and A549 lysates using STAT1 antibody.



Immunoprecipitation analysis of STAT1 (6C2) in HeLa lysates using STAT1 antibody.