


KDM5C (2E4) Mouse mAb

ABZ3972


Product Overview

Name	KDM5C (2E4) Mouse mAb
Accession(Primary)	P41229
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 176 kDa; Observed MW: 220 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (P41229)

Gene Name KDM5C**About Gene Symbol**

Western blot analysis of KDM5C / Jarid1C / SMCX in HeLa and 293 lysates using KDM5C / Jarid1C / SMCX antibody.



Immunocytochemistry analysis of KDM5C in HeLa cells using KDM5C antibody.



Western blot analysis of KDM5C (2E4) in T47D, MCF-7 and HeLa lysates using KDM5C antibody.