

STAT3 (6H10) Mouse mAb


ABZ3973

Product Overview


Name	STAT3 (6H10) Mouse mAb
Accession(Primary)	P40763
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 79,86 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P40763)


Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of STAT3 (6H10) in HeLa, 3T3, PC-12 lysates using STAT3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma using STAT3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of STAT3 (6H10) in Jurkat, Raw264.7, Ramos, A549 lysates using Stat3 (6H1) antibody.

