

STAT3 (4C8) Mouse mAb


ABZ3974

Product Overview

Name	STAT3 (4C8) Mouse mAb
Accession(Primary)	P40763
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat, Hamster
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P40763)

Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of STAT3 in
CHOK1,MCF7,A549,C6,3T3,K562 and HeLa lysates
using STAT3 antibody.