


MLH1 (4C1) Mouse mAb

ABZ3975


Product Overview

Name	MLH1 (4C1) Mouse mAb
Accession(Primary)	P40692
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (P40692)

Gene Name MLH1**About Gene Symbol**

Immunocytochemistry analysis of MLH1 (4C1) in HeLa cells using MLH1 antibody.



Western blot analysis of MLH1 (4C1) in HeLa, Jurkat, 293T and SW480 lysates using MLH1 (4C1) antibody



Western blot analysis of MLH1 (4C1) in SW48, 293T, Jurkat and HeLa lysates using MLH1 (4C1) antibody