

FEN1 (1E7) Mouse mAb


ABZ3977

Product Overview


Name	FEN1 (1E7) Mouse mAb
Accession(Primary)	P39748
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P39748)


Gene Name	FEN1
About Gene Symbol	FEN1 (Flap Structure-Specific Endonuclease 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of FEN1 in HeLa, Jurkat, 3T3, C6, CEM and Ramos lysates using FEN1 antibody.



Immunoprecipitation analysis of FEN1 (1E7) in HeLa lysates using FEN1 antibody.



Western blot analysis of FEN1 (1E7) in HeLa, Jurkat, K562 and 3T3 lysates using FEN1 antibody.