


PP1C gamma (2F2A3) Mouse mAb

ABZ3978

Product Overview

Name	PP1C gamma (2F2A3) Mouse mAb
Accession(Primary)	P36873
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P36873)

Gene Name PPP1CC**About Gene Symbol**

Western blot analysis of PP1C in HeLa, 293T and Jurkat lysates using PP1C antibody.

