

beta Catenin Mouse mAb


ABZ3980

Product Overview


Name	beta Catenin Mouse mAb
Accession(Primary)	P35222
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Zebrafish
Calculated MW	Calculated MW: 85 kDa; Observed MW: 92 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC-1:200-1:500; IF-1:50-1:200

Target information (P35222)


Gene Name	CTNNB1
About Gene Symbol	



Western blot analysis of beta Catenin (7E2) in HeLa, 293T, mouse Liver, rat Liver lysates using β Catenin antibody.



Western blot analysis of beta Catenin (7E2) in T47D, 3T3, COS7 and C6 lysates using β Catenin antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using beta Catenin (7E2) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.