


beta Catenin Mouse mAb

ABZ3981

Product Overview

Name	beta Catenin Mouse mAb
Accession(Primary)	P35222
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 92 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P35222)

Gene Name CTNNB1**About Gene Symbol**

Western blot analysis of β Catenin in HCT116, C6 and 3T3 lysates using β Catenin antibody.

