

MPI (4G9) Mouse mAb

ABZ3982


Product Overview

Name	MPI (4G9) Mouse mAb
Accession(Primary)	P34949
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (P34949)

Gene Name MPI

About Gene Symbol



Western blot detection of Mannose Phosphate Isomerase in rat kidney, rat brain, A549 and Lncap lysates using Mannose Phosphate Isomerase mouse mAb (1:1000 diluted).



Immunocytochemistry stain of HeLa using
Mannose Phosphate Isomerase mouse mAb
(1:300).