

CD68 (9H5) Mouse mAb

ABZ3983


Product Overview

Name	CD68 (9H5) Mouse mAb
Accession(Primary)	P34810
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:200-1:400; IF-1:50-1:200


Target information (P34810)

Gene Name CD68


About Gene Symbol




Immunohistochemistry analysis of paraffin-embedded Human Tonsil tissue using CD68 (9H5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.




Immunohistochemical analysis of paraffin-embedded Human tonsils using CD68 (9H5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunofluorescence analysis of CD68 (9H5) in Human lung cancer tissue using CD68 (9H5) antibody (red) and DAPI (blue).



Immunofluorescence analysis of CD68 (9H5) in mouse spleen tissue using CD68 (9H5) antibody (6F3) (red), and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using CD68 (9H5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.