


Histone H2B (7E2) Mouse mAb

ABZ3984


Product Overview

Name	Histone H2B (7E2) Mouse mAb
Accession(Primary)	P33778
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P33778)

Gene Name H2BC3**About Gene Symbol**

Western blot analysis of Histone H2B in HeLa and COS7 lysates using Histone H2B antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma using Histone H2B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.