

MCM3 Rabbit mAb


ABZ3987

Product Overview


Name	MCM3 Rabbit mAb
Accession(Primary)	P25205
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 91 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; IHC-P-1:50-1:100

Target information (P25205)


Gene Name	MCM3
About Gene Symbol	MCM3 (Minichromosome Maintenance Complex Component 3) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic



Western blot analysis of MCM3 in Jurkat, C6, HeLa lysates using MCM3 antibody.



Immunocytochemistry analysis of MCM3(green) in HeLa using MCM3 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MCM3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.