

Cyclin D2 (6E11) Mouse mAb

ABZ3988


Product Overview

Name	Cyclin D2 (6E11) Mouse mAb
Accession(Primary)	P30279
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 33 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P30279)

Gene Name CCND2

About Gene Symbol



Western blot analysis of Cyclin D2 in CHO-K1 lysates and CHO-K1 transfected by Cyclin D2 fragment EGFP fusion protein lysates using Cyclin D2 antibody.

