

## eNOS (2E8) Mouse mAb

ABZ3989


### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | eNOS (2E8) Mouse mAb                         |
| <b>Accession(Primary)</b>   | P29474                                       |
| <b>Host Species</b>         | Mouse                                        |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG1                                         |
| <b>Species Reactivity</b>   | Human, Mouse, Rat, Rabbit                    |
| <b>Calculated MW</b>        | Calculated MW: 133 kDa; Observed MW: 133 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:2000                              |

### Target information (P29474)

**Gene Name** NOS3

**About Gene Symbol**



Western blot analysis of eNOS (2E8) in rat Heart, huvec , Jurkat lysates using eNOS (2E8) antibody.

