


SHP1 (2B7) Mouse mAb

ABZ3990


Product Overview

Name	SHP1 (2B7) Mouse mAb
Accession(Primary)	P29350
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P29350)

Gene Name PTPN6**About Gene Symbol**

Western blot analysis of SHP1 in Human tonsil, Daudi, Jurkat and Raji lysates using SHP1 antibody.



Immunoprecipitation analysis of SHP1 (2B7) in
Raji lysates using SHP1 antibody.