

ERK1/2 (4A4) Mouse mAb


ABZ3993

Product Overview

Name	ERK1/2 (4A4) Mouse mAb
Accession(Primary)	P27361
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44,42 kDa; Observed MW: 42,44 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P27361)

Gene Name	MAPK3/MAPK1
About Gene Symbol	MAPK3 (Mitogen-Activated Protein Kinase 3) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ERK1/2 (4A4) in EGFtreated A549 , PDGFtreated 3T3 and TPAtreated 3T3 lysates using ERK1/2 (4A4) antibody.