

RPA70 (8C3) Mouse mAb


ABZ3997

Product Overview


Name	RPA70 (8C3) Mouse mAb
Accession(Primary)	P27694
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P27694)


Gene Name	RPA1
About Gene Symbol	RPA1 (Replication Protein A1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of RPA70 in HeLa, A431, MCF-7, COS7, HT-1080 and K562 lysates using RPA70 antibody.




Immunocytochemistry analysis of RPA70 (8C3) in HeLa using RPA70 antibody.



Immunoprecipitation analysis of RPA70 (8C3) in HeLa lysates using RPA7 antibody.



Western blot analysis of RPA70 (8C3) in HeLa, 293T, C6, 3T3 and K562 lysates using RPA70 (8C3) antibody



Western blot analysis of RPA70 in HUVEC, Hacat, 3T3, MCF-7 and A549 lysates using RPA70 antibody.