

RPA70 Mouse mAb


ABZ3998

Product Overview

Name	RPA70 Mouse mAb
Accession(Primary)	P27694
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P27694)

Gene Name	RPA1
About Gene Symbol	RPA1 (Replication Protein A1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of RPA70 in C6, COS7, K562 and HeLa lysates using RPA70 antibody.