


S100A4 (3D10) Mouse mAb

ABZ4001


Product Overview

Name	S100A4 (3D10) Mouse mAb
Accession(Primary)	P26447
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P26447)

Gene Name S100A4**About Gene Symbol**

Western blot analysis of S100A4 in HeLa and A549 lysates using S100A4 antibody.



Immunoprecipitation analysis of S100A4 (3D10) in HeLa lysates using S100A4 antibody.