

MCM3 Rabbit mAb

ABZ4003

Product Overview

Name	MCM3 Rabbit mAb
Accession(Primary)	P25205
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	IHC-1:100-1:200

Target information (P25205)

Gene Name	MCM3
About Gene Symbol	MCM3 (Minichromosome Maintenance Complex Component 3) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic