

Desmin (4H1) Mouse mAb


ABZ4007

Product Overview


Name	Desmin (4H1) Mouse mAb
Accession(Primary)	P17661
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (P17661)

Gene Name	DES
About Gene Symbol	DES (Desmin) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cardiovascular disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded mouse Skeletal Muscle Tissue using Desmin (4H1) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Desmin (4H1) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.