

EpCAM (1D5) Mouse mAb

ABZ4008


Product Overview

Name	EpCAM (1D5) Mouse mAb
Accession(Primary)	P16422
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 35 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20


Target information (P16422)

Gene Name EPCAM

About Gene Symbol



Western blot analysis of EpCAM (1D5) in CHO-K1 transfected by FlagEpCAM, MCF-7(EpCAM positive), HCT116(EpCAM positive), and HeLa(EpCAM negative) lysates using EpCAM (1D5) antibody.



Immunoprecipitation analysis of EpCAM (1D5) in HCT116 and HeLa lysates using EpCAM antibody.