


Hsp27 (7E5) Mouse mAb

ABZ4009


Product Overview

Name	Hsp27 (7E5) Mouse mAb
Accession(Primary)	P04792
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P04792)

Gene Name HSPB1**About Gene Symbol**

Western blot analysis of Hsp27 in HeLa, A431, COS7, Lncap and MCF-7 lysates using Hsp27 antibody.



Immunocytochemistry analysis of Hsp27 (7E5) in Hela using Hsp27 antibody.