


Hsp27 (9G1) Mouse mAb

ABZ4010


Product Overview

Name	Hsp27 (9G1) Mouse mAb
Accession(Primary)	P04792
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 23 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P04792)

Gene Name HSPB1**About Gene Symbol**

Western blot analysis of Hsp27 (9G1) in HeLa lysates using Hsp27 (9G1) antibody.



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma using HSP27 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Hsp27 (9G1) in K562, HeLa lysates using HSP27 (9G1) antibody.