


S100 beta (5G3) Mouse mAb

ABZ4013

Product Overview

Name	S100 beta (5G3) Mouse mAb
Accession(Primary)	P04271
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P04271)

Gene Name S100B**About Gene Symbol**

Western blot analysis of S100B (5G3) in rat Brain lysates using S100B antibody.

