

SOD2 (4C4) Mouse mAb


ABZ4014

Product Overview

Name	SOD2 (4C4) Mouse mAb
Accession(Primary)	P04179
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 25 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P04179)

Gene Name	SOD2
About Gene Symbol	SOD2 (Superoxide Dismutase 2) is a metabolic enzyme involved in cellular metabolism. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SOD2 in HeLa, THP-1, A549 Ramos and K562 lysates using SOD2 antibody.