

ER alpha Mouse mAb

ABZ4015


Product Overview

Name	ER alpha Mouse mAb
Accession(Primary)	P03372
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 66 kDa; Observed MW: 66 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P03372)

Gene Name ESR1

About Gene Symbol



Western blot analysis of Estrogen Receptor alpha (6F11) in CHO-K1 (B) and CHO-K1 transfected by Estrogen Receptor α fragment(A) lysates using Estrogen Receptor α antibody.

