


Transferrin Mouse mAb

ABZ4016

Product Overview

Name	Transferrin Mouse mAb
Accession(Primary)	P02787
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; ELISA-1:10000

Target information (P02787)

Gene Name TF**About Gene Symbol**

Western blot analysis of Transferrin (1G7) in Human Serumlysates using Transferrin (1G7) antibody.

