

Human Serum Albumin (1A8) Mouse mAb


ABZ4018

Product Overview


Name	Human Serum Albumin (1A8) Mouse mAb
Accession(Primary)	P02768
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 69 kDa; Observed MW: 69 kDa
Recommended Dilution	WB-1:500-1:1000; ELISA-1:10000

Target information (P02768)


Gene Name	ALB
About Gene Symbol	ALB (Albumin) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Human Serum Albumin in 0.5nl Human serum and 5ng ALB Recombinant antigens lysates using Human Serum Albumin antibody.



Indirect ELISA assay for antiHuman Serum Albumin antibody. Antigen coating concentration: 8ug/ml.



Standard Curve for Human Serum Albumin (Analyte: Human Serum Albumin protein); using Capture antibody mouse monoclonal [1A8C5B4] to Human Serum Albumin at 8ug/ml and Detector antibody Rabbit polyclonal antiserum to Human Serum Albumin.