

FN Mouse mAb


ABZ4021

Product Overview


Name	FN Mouse mAb
Accession(Primary)	P02751
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 272 kDa; Observed MW: 285 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:100-1:500; IHC-1:50-1:100

Target information (P02751)

Gene Name	FN1
About Gene Symbol	FN1 (Fibronectin 1) is a protein kinase and signaling regulator involved in cell migration and invasion. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Fibronectin (1E2) in MCF-7 and Hela lysates using Fibronectin antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney using Fibronectin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.