

C Reactive Protein (1H3) Mouse mAb

ABZ4022


Product Overview

Name	C Reactive Protein (1H3) Mouse mAb
Accession(Primary)	P02741
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000


Target information (P02741)

Gene Name CRP

About Gene Symbol



sELISA standard Curve for CReactive Protein:
 Capture Antibody mouse mAb to CReactive Protein
 at 1ug/ml and 2012781H3 was used for detecting.



sELISA standard Curve for CReactive Protein:
 Capture Antibody Rabbit pAb to CReactive Protein
 at 4ug/ml and 2012781H3 was used for detecting.