

C Reactive Protein (6D1) Mouse mAb

ABZ4023


Product Overview

Name	C Reactive Protein (6D1) Mouse mAb
Accession(Primary)	P02741
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000

Target information (P02741)

Gene Name CRP

About Gene Symbol



sELISA standard Curve for CReactive Protein:
 Capture Antibody mouse mAb(2012786D1) to
 CReactive Protein at 2ug/ml and 2012781H3 was
 used for detecting.