


Apolipoprotein A I Mouse mAb

ABZ4024

Product Overview

Name	Apolipoprotein A I Mouse mAb
Accession(Primary)	P02647
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 31 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P02647)

Gene Name APOA1**About Gene Symbol**

Western blot analysis of ApoA1 (2D6) in HepG2 and Human serum lysates using ApoA1 (2D6) antibody.

