

## Collagen IV alpha 1 (2G3) Mouse mAb

ABZ4025


### Product Overview

|                             |                                     |
|-----------------------------|-------------------------------------|
| <b>Name</b>                 | Collagen IV alpha 1 (2G3) Mouse mAb |
| <b>Accession(Primary)</b>   | P02462                              |
| <b>Host Species</b>         | Mouse                               |
| <b>Clonality</b>            | Monoclonal                          |
| <b>Isotype</b>              | IgG1                                |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                   |
| <b>Recommended Dilution</b> | IHC-1:50-1:100                      |


### Target information (P02462)

**Gene Name** COL4A1

**About Gene Symbol**



Immunohistochemistry analysis of paraffin-embedded Human Liver Carcinoma Tissue using Collagen IV alpha 1 (2G3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Collagen IV alpha 1 (2G3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.