

Collagen I alpha 1 Mouse mAb

ABZ4026


Product Overview

Name	Collagen I alpha 1 Mouse mAb
Accession(Primary)	P02452
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100


Target information (P02452)

Gene Name COL1A1

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma Tissue using Collagen I alpha 1 (4C4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Collagen I alpha 1 (4C4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.