


Human IgM (3E12) Mouse mAb

ABZ4027


Product Overview

Name	Human IgM (3E12) Mouse mAb
Accession(Primary)	P01871
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Recommended Dilution	IHC-1:50-1:100; ELISA-1:10000

Target information (P01871)

Gene Name IGHM**About Gene Symbol**

Indirect ELISA assay for mouse Human IgM antibody. Antigen coating concentration: 2ug/ml.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Human IgM (3E12) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.