

TNF alpha (9D9) Mouse mAb

ABZ4031

Product Overview


Name	TNF alpha (9D9) Mouse mAb
Accession(Primary)	P01375
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 17,25 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P01375)


Gene Name TNF

About Gene Symbol

Recombinant Human TNF a Protein



Western blot analysis of TNF alpha in Recombinant Human TNF a Protein using TNF alpha antibody.



Immunohistochemistry analysis of paraffin-embedded Human Tonsil Tissue using TNF alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.