


c-Fos (7D6) Mouse mAb

ABZ4036

Product Overview

Name	c-Fos (7D6) Mouse mAb
Accession(Primary)	P01100
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Rat, Cat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 62 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P01100)

Gene Name FOS**About Gene Symbol**

Western blot analysis of c-Fos in HCT116, COS7 and C6 lysates using c-Fos antibody.

