

Cystatin C (3B11) Mouse mAb


ABZ4038

Product Overview

Name	Cystatin C (3B11) Mouse mAb
Accession(Primary)	P01034
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	ELISA-1:10000

Target information (P01034)

Gene Name	CST3
About Gene Symbol	CST3 (Cystatin C) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Standard Curve for Human Cystatin C: Capture Antibody mouse mAb [6F12C7D8] to Human Cystatin C at 4ug/ml and Detector Antibody mouse mAb [2009125C3:5C3C7.2009123B11:3B11E9] to Human Cystatin C at 0.5ug/ml.