

EGFR (5E9) Mouse mAb


ABZ4042

Product Overview


Name	EGFR (5E9) Mouse mAb
Accession(Primary)	P00533
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 134 kDa; Observed MW: 170 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P00533)


Gene Name	EGFR
About Gene Symbol	EGFR (Epidermal Growth Factor Receptor) is a receptor tyrosine kinase involved in immune signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of EGFR (5E9) in A431, HeLa lysates using EGFR (5E9) antibody.



Immunoprecipitation analysis of EGFR (5E9) in HeLa lysates using EGFR (5E9) antibody



Western blot analysis of EGFR (5E9) in EGF
Receptor/ErbB1 in HeLa lysates using EGF
Receptor/ErbB1 (5E9) antibody.