


JAK2 Mouse mAb

ABZ4050


Product Overview

Name	JAK2 Mouse mAb
Accession(Primary)	O60674
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (O60674)

Gene Name JAK2**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK2 (7H5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using JAK2 (7H5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.