


KDM1/LSD1 Mouse mAb

ABZ4051

Product Overview

Name	KDM1/LSD1 Mouse mAb
Accession(Primary)	O60341
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 93 kDa; Observed MW: 93 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O60341)

Gene Name KDM1A**About Gene Symbol**

Western blot analysis of KDM1 / LSD1 in A549 and Jurkat lysates using KDM1 / LSD1 antibody.

