


KDM1/LSD1 Mouse mAb

ABZ4052


Product Overview

Name	KDM1/LSD1 Mouse mAb
Accession(Primary)	O60341
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 93 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20


Target information (O60341)

Gene Name KDM1A**About Gene Symbol**

Western blot analysis of KDM1/LSD1 in HeLa, HCT116 and COS7 lysates using KDM1/LSD1 antibody.



Immunocytochemistry analysis of KDM1A (1E5) in HeLa using KDM1/LSD1 antibody.



Immunoprecipitation analysis of KDM1A (1E5) in HeLa lysates using KDM1/LSD1 antibody.