


## CD133 (6H10) Mouse mAb


ABZ4053

### Product Overview

|                             |                                             |
|-----------------------------|---------------------------------------------|
| <b>Name</b>                 | CD133 (6H10) Mouse mAb                      |
| <b>Accession(Primary)</b>   | O43490                                      |
| <b>Host Species</b>         | Mouse                                       |
| <b>Clonality</b>            | Monoclonal                                  |
| <b>Isotype</b>              | IgG2a                                       |
| <b>Species Reactivity</b>   | Human                                       |
| <b>Calculated MW</b>        | Calculated MW: 97 kDa; Observed MW: 133 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                             |

### Target information (O43490)

**Gene Name** PROM1**About Gene Symbol**



Western blot analysis of CD133 (6H10) in CHO-K1 cells nontransfected (A) or transfected (B) with CD133 using CD133 antibody.