

TAK1 Mouse mAb


ABZ4055

Product Overview

Name	TAK1 Mouse mAb
Accession(Primary)	O43318
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 67 kDa; Observed MW: 67,78 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O43318)

Gene Name	MAP3K7
About Gene Symbol	MAP3K7 (Mitogen-Activated Protein Kinase Kinase Kinase 7) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and TGF-β signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TAK1 in K562, HeLa, MCF-7, COS7, 3T3 and C6 lysates using TAK1 antibody.