


FGFR4 Mouse mAb

ABZ4057


Product Overview

Name	FGFR4 Mouse mAb
Accession(Primary)	P22455
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P22455)

Gene Name FGFR4**About Gene Symbol**

Western blot analysis of FGFR4 in Human FGFR4 (aa22369) lysates using FGFR4 antibody.



Immunofluorescence analysis of FGFR4 in HEK293 cells trasfected with FGFR4hlgGFc using FGFR4 antibody(green) Blue: DRAQ5 fluorescent DNA dye.