

CD23 (3D1) Mouse mAb

ABZ4059


Product Overview

Name	CD23 (3D1) Mouse mAb
Accession(Primary)	P06734
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100; IF-1:50-1:200


Target information (P06734)

Gene Name FCER2


About Gene Symbol




Immunofluorescence analysis of CD23 (3D1) in Human stomach using CD23 (3D1) antibody (red), and DAPI (blue).




Immunofluorescence analysis of CD23 (3D1) in rat lung tissue using CD23 antibody(1E9)(red),and DAPI (blue).




Immunohistochemistry analysis of paraffin-embedded Human lung tissue using CD23 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded rat kidney tissue using CD23 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded mouse lung tissue using CD23 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using CD23 (3D1) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.