


CD15 Mouse mAb

ABZ4062


Product Overview

Name	CD15 Mouse mAb
Accession(Primary)	P22083
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	IHC-1:50-1:100; IF-1:50-1:200


Target information (P22083)

Gene Name FUT4**About Gene Symbol**


Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using CD15 (8C10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.




Immunofluorescence analysis of CD15 (8C10) in Human livercancer tissue using CD15 antibody (red), and DAPI (blue).




Immunohistochemistry analysis of paraffin-embedded rat Brain Tissue using CD 15 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human stomach using CD15 (8C10) antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human stomach using CD15 (8C10) antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human stomach using CD15 (8C10) antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.