

ITGA2B/CD41 (9G4) Mouse mAb

ABZ4063


Product Overview

Name	ITGA2B/CD41 (9G4) Mouse mAb
Accession(Primary)	P08514
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	IHC-1:50-1:100


Target information (P08514)

Gene Name ITGA2B

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human breastcancer tissue using Integrin alpha lib (9G4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Integrin alpha lib (9G4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.